

SARS-CoV-2/Flu A/Flu B Ag Rapid Test



SARS-CoV-2/Flu A/Flu B Ag Rapid Test

VivaDiag™ SARS-CoV-2/Flu A/Flu B Ag Rapid Test is for the simultaneous detection and differentiation of the antigens from SARS-CoV-2, influenza A and influenza B in the specimens from individual suspected of respiratory viral infections. The test will serve as a powerful weapon for the healthcare professionals to fight against the COVID-19 pandemic.


Specification

Test Parameters	SARS-CoV-2 Ag /Flu A Ag /Flu B Ag
Test Principle	Immunochromatography
Sample Type	Nasal swab / Nasopharyngeal swab / Oropharyngeal swab
Sample Volume	3 drops
Test Time	15 min
Operation Temperature	15-30°C
Storage Temperature	2-30°C
Shelf Life	24 months




Test Procedure


Specimen Collection



Solution Preparation




Detection



3. Read the test result

Result Interpretation



Note: The VivaDiag™ SARS-CoV-2/Flu A/Flu B Ag Rapid Test has ONLY been designed to function as a supplementary tool in conjunction with nucleic acid test for the diagnosis of suspected cases of SARS-CoV-2, influenza A or influenza B. Thus, the results from rapid test should not be used as the basis to diagnose or exclude SARS-CoV-2, influenza A or influenza B infection. Furthermore, the test results should not be used to determine or inform the infection status.

Plusmedics Supplies Sdn Bhd

Tel: +604 - 4121 292 / +6012-5661230

Email: kiru@plusmedics.health

Website: www.plusmedics.health

Medical Device
AUTHORITY

MALAYSIA